



HEART RATE ESTIMATOR



ECG Signal



```
clc;
clear;
close all;

% Sampling frequency
fs = 500; % Standard ECG sampling rate (Hz)
t_total = 5; % Total duration (seconds) for 5 beats
t = 0:1/fs:t_total; % Time vector

% Heart rate (bpm)
HR = 60; % 60 BPM (1 beat per second)
f_HR = HR / 60; % Convert BPM to Hz
ecg_period = 1 / f_HR; % One heartbeat duration

% ECG waveform parameters (P-QRS-T complex)
A_p = 0.1; % P wave amplitude
A_q = -0.15; % Q wave amplitude
A_r = 1.2; % R wave amplitude (main peak)
A_r_neg = -0.5; % Occasional Negative R peak
A_s = -0.4; % S wave amplitude
A_t = 0.3; % T wave amplitude


% Wave timing (based on ECG characteristics)
T_p = -0.2 * ecg_period; % P wave occurs before QRS
T_q = -0.05 * ecg_period; % Q wave slightly before R peak
T_r = 0; % R wave at center
T_r_neg = 0.02 * ecg_period; % Negative R peak (occasional)
T_s = 0.05 * ecg_period; % S wave after R
T_t = 0.3 * ecg_period; % T wave occurs later
```

```


% Wave durations (Gaussian standard deviations)
sigma_p = 0.08 * ecg_period; % P wave width
sigma_q = 0.02 * ecg_period; % Q wave width
sigma_r = 0.03 * ecg_period; % R wave width
sigma_r_neg = 0.03 * ecg_period; % Negative R wave width
sigma_s = 0.02 * ecg_period; % S wave width
sigma_t = 0.12 * ecg_period; % T wave width
% Generate ECG for one heartbeat
t_single = -ecg_period/2:1/fs:ecg_period/2;
ecg_heartbeat = zeros(size(t_single));
for i = 1:length(t_single)
    ecg_heartbeat(i) = ...
        A_p * exp(-((t_single(i) - T_p).^2) / (2 * sigma_p^2)) + ... % P wave
        A_q * exp(-((t_single(i) - T_q).^2) / (2 * sigma_q^2)) + ... % Q wave
        A_r * exp(-((t_single(i) - T_r).^2) / (2 * sigma_r^2)) + ... % R wave
        (mod(i, round(fs * ecg_period * 2)) == 0) * A_r_neg * exp(-((t_single(i) - T_r_neg).^2) / (2 * sigma_r_neg^2)) + ... % Occasional Negative R peak
        A_s * exp(-((t_single(i) - T_s).^2) / (2 * sigma_s^2)) + ... % S wave
        A_t * exp(-((t_single(i) - T_t).^2) / (2 * sigma_t^2)); % T wave
end
% Generate a full ECG signal with 5 beats
num_beats = 5; % Exactly 5 R-peaks
final_ecg = [];
final_t = [];
for i = 1:num_beats
    final_ecg = [final_ecg, ecg_heartbeat]; % Append each heartbeat
    final_t = [final_t, t_single + (i-1) * ecg_period]; % Adjust time axis
end
% Add Moderate Noise Components
baseline_wander = 0.05 * sin(2 * pi * 0.5 * final_t); % Mild baseline drift (0.5Hz)
muscle_noise = 0.02 * randn(size(final_t)); % Smaller high-frequency noise
powerline_noise = 0.01 * sin(2 * pi * 50 * final_t); % Weaker 50Hz powerline interference
% Final Noisy ECG Signal
ecg_noisy = final_ecg + baseline_wander + muscle_noise + powerline_noise;
% Save ECG data for further processing
fileID = fopen('ecg_noisy.txt', 'w');
fprintf(fileID, '%f\n', ecg_noisy);
fclose(fileID);
plot(final_t, ecg_noisy, 'r', 'LineWidth', 1.5);
title('ECG Signal');
xlabel('Time (s)');
ylabel('Amplitude (mV)');
grid on;
axis tight;

```


Filter




Differentiator



Rectifier



Integrator



Verilog Code

```
module heart_rate_pipeline (  
    input clk,  
    input rst_n,  
    input signed [7:0] Xin,  
    output signed [12:0] Yout  
);  
    // Stage 1: Differentiation (Xin - Xin_prev)  
    reg signed [7:0] Xin_Reg [1:0];  
    always @(posedge clk or negedge rst_n) begin  
        if (!rst_n) begin  
            Xin_Reg[0] <= 8'd0;  
            Xin_Reg[1] <= 8'd0;  
        end else begin  
            Xin_Reg[1] <= Xin_Reg[0];  
            Xin_Reg[0] <= Xin;  
        end  
    end  
    wire signed [8:0] diff_out = Xin_Reg[0] - Xin_Reg[1];  
  
    // Stage 2: Rectification (Absolute Value)  
    reg signed [8:0] rect_out;  
    always @(posedge clk or negedge rst_n) begin  
        if (!rst_n)  
            rect_out <= 9'd0;  
        else  
            rect_out <= (diff_out[8] == 1'b0) ? diff_out : (~diff_out + 1'b1);  
        end  
    end
```

```

// Stage 3: Integration (16-element shift register + 4-level adder tree)
reg signed [8:0] Xin_Reg_Inte [15:0];
integer i; // Declare integer outside the always block
always @(posedge clk or negedge rst_n) begin
    if (!rst_n) begin
        for (i = 0; i < 16; i = i + 1)
            Xin_Reg_Inte[i] <= 9'd0;
    end else begin
        for (i = 15; i > 0; i = i - 1)
            Xin_Reg_Inte[i] <= Xin_Reg_Inte[i - 1];
        Xin_Reg_Inte[0] <= rect_out;
    end
end

// 4-level adder tree for accumulation
wire signed [12:0] sum_lvl_2;
assign sum_lvl_2 = Xin_Reg_Inte[0] + Xin_Reg_Inte[1] + Xin_Reg_Inte[2] +
Xin_Reg_Inte[3] +
                Xin_Reg_Inte[4] + Xin_Reg_Inte[5] + Xin_Reg_Inte[6] + Xin_Reg_Inte[7] +
                Xin_Reg_Inte[8] + Xin_Reg_Inte[9] + Xin_Reg_Inte[10] + Xin_Reg_Inte[11] +
                Xin_Reg_Inte[12] + Xin_Reg_Inte[13] + Xin_Reg_Inte[14] + Xin_Reg_Inte[15];
assign Yout = sum_lvl_2;
endmodule

```

Testbench for Filter :

```

`timescale 1ns/1ps
module heart_rate_pipeline_tb;
    // Inputs
    reg clk;
    reg rst_n;
    reg signed [7:0] Xin;
    // Outputs
    wire signed [12:0] Yout;
    // Instantiate the Unit Under Test (UUT)
    heart_rate_pipeline uut (
        .clk(clk),
        .rst_n(rst_n),
        .Xin(Xin),
        .Yout(Yout)
    );

```

```

initial begin
    clk = 0;
    forever #5 clk = ~clk; // 10ns period clock
end
initial begin
    rst_n = 0;
    #20; // Hold reset for 20ns
    rst_n = 1;
end
// Read ECG data from file and apply to Xin
integer file;
integer scan_result;
real ecg_value;
initial begin
    file = $fopen("ecg_signal.txt", "r");
    if (file == 0) begin
        $display("Error: Could not open file.");
        $finish;
    end
// Initialize Xin
    Xin = 0;
    #30;
    while (!$feof(file)) begin
        scan_result = $fscanf(file, "%f\n", ecg_value); // Read a floating-point value from the file
        if (scan_result != 1) begin
            $display("Error: Failed to read ECG value from file.");
            $finish;
        end
        Xin = $rtoi(ecg_value * 127); // Convert to 8-bit fixed-point
        #10;
    end
//Close the file
    $fclose(file);


    // Finish the simulation
    #100;
    $finish;
end

// Monitor the outputs
initial begin
    $monitor("Time: %0t | Xin: %d | Yout: %d", $time, Xin, Yout);
    $dumpfile("heart_rate_pipeline");
    $dumpvars(0, heart_rate_pipeline_tb);
end

endmodule

```


Filtered Output :



Peak Detector

Features

**Threshold
Adjustment for
Peak detection**



**Fixed RR interval
to avoid Fault
Detection**

Code :

```

module self_adaptive_threshold (
    input clk,
    input rst_n,
    input signed [12:0] transformed_signal,
    output reg peak_detected);
    reg signed [12:0] signal_buffer [0:255];
    reg [7:0] index;
    integer j;
    always @(posedge clk or negedge rst_n) begin
        if (!rst_n)
            index <= 8'd0;
        else if (index < 8'd255)
            index <= index + 8'd1;
        else
            index <= 8'd0;
    end
    always @(posedge clk or negedge rst_n) begin
        if (!rst_n) begin
            for (j = 0; j < 256; j = j + 1)
                signal_buffer[j] <= 13'd0;
        end else
            signal_buffer[index] <= transformed_signal;
    end
end


```

```

// Calculate the maximum value in the buffer
reg signed [12:0] max_value, temp_max_value;
always @(posedge clk or negedge rst_n) begin
    if (!rst_n)
        max_value <= 13'd0;
    else begin
        temp_max_value = 13'd0;
        for (j = 0; j < 256; j = j + 1)
            if (signal_buffer[j] > temp_max_value)
                temp_max_value = signal_buffer[j];
        max_value <= temp_max_value;
    end
end
// Set the dynamic threshold
reg signed [12:0] threshold;
always @(posedge clk or negedge rst_n) begin
    if (!rst_n)
        threshold <= 13'd0;
    else
        threshold <= (max_value >> 1) + 13'd50; // Use arithmetic right shift
end
// Detect peaks
always @(posedge clk or negedge rst_n) begin
    if (!rst_n)
        peak_detected <= 1'b0;
    else
        peak_detected <= (transformed_signal > threshold);
end
endmodule

```

Output :



RR Interval Calculator :

Average of N intervals is considered as the RR interval of the ECG to enhance accuracy

Code :

```
module RR_Interval_Calculator (
    input clk,
    input rst_n,
    input peak_detected,
    output reg [15:0] avg_interval
);
    parameter NUM_PEAKS = 4;
    reg output_valid;
    reg [15:0] current_time;
    reg [15:0] sum_intervals;
    reg [3:0] index;
    reg [15:0] current_interval;
    reg [15:0] prev_time;
    reg calibration_done;
    reg prev_peak_detected;
    wire peak_rising_edge;
    assign peak_rising_edge = peak_detected && !prev_peak_detected;
    integer i;
    always @(posedge clk or negedge rst_n) begin
        if (!rst_n) begin
            prev_time = 16'd0;
            current_time = 16'd0;
            index = 4'd0;
            sum_intervals = 16'd0;
            avg_interval = 16'd0;
            calibration_done = 1'b0;
            output_valid = 1'b0;
            prev_peak_detected = 1'b0;
        end
    end
end
```

```

else begin
    prev_peak_detected = peak_detected;
    if (peak_rising_edge) begin
        if (prev_time == 16'd0) begin
            prev_time = current_time;
        end
    else begin
        current_interval = current_time - prev_time;
        sum_intervals = sum_intervals + (current_time - prev_time);
        prev_time = current_time;
        if (!calibration_done) begin
            if (index == NUM_PEAKS-1) begin
                avg_interval = sum_intervals / (NUM_PEAKS);
                calibration_done = 1'b1;
                output_valid = 1'b1;
            end else begin
                index = index + 1;
            end
        end
    end
end
end
else begin
    current_time = current_time + 1;
end
end
end
endmodule

```

Testbench :

```

`timescale 1ns/1ps
module RR_interval_calculator_tb;
    reg clk;
    reg rst_n;
    reg signed [7:0] Xin;
    wire signed [12:0] Yout;
    wire peak_detected;
    wire [15:0] avg_RR_interval;

```

```

heart_rate_pipeline uut_pipeline (
    .clk(clk),
    .rst_n(rst_n),
    .Xin(Xin),
    .Yout(Yout)
);
self_adaptive_threshold uut_threshold (
    .clk(clk),
    .rst_n(rst_n),
    .transformed_signal(Yout),
    .peak_detected(peak_detected)
);
RR_Interval_Calculator uut_rr_calculator (
    .clk(clk),
    .rst_n(rst_n),
    .peak_detected(peak_detected),
    .avg_interval(avg_RR_interval)
);
initial begin
    clk = 0;
end
always #10 clk = ~clk;
initial begin
    rst_n = 0;
    #20;
    rst_n = 1;
end
integer file;
integer scan_result;
real ecg_value;
initial begin
    file = $fopen("ecg_signal.txt", "r");
    if (file == 0) begin
        $display("Error: Could not open file.");
        $finish;
    end
    Xin = 0;
    #30;
    while (!$feof(file)) begin
        scan_result = $fscanf(file, "%f\n", ecg_value);
        if (scan_result != 1) begin
            $display("Error: Failed to read ECG value from file.");
            $finish;
        end
        Xin = $rtoi(ecg_value * 128);
        #10;
    end
    $fclose(file);
    #100;
    $finish;
end

```

```

initial begin
    $monitor("Time: %0t | Xin: %d | Yout: %d | Peak Detected: %b | Avg RR Interval: %d",
    $time, Xin, Yout, peak_detected, avg_RR_interval);
    $dumpfile("RR_interval_calculator.vcd");
    $dumpvars(0, RR_interval_calculator_tb);
end

endmodule

```

Output :

